

## Transfer Abutment and Transfer Cap

### OLD

Transfer section		Cone
Gingiva / Emergence profile		Round
Connection		HEX

Can no longer be ordered

### NEW

	Cone with locking function
	Round with 2 faces
	Prism-HEX

	GH 4.0 mm	GH 7.0 mm
3.3 mm	936024	936025

### OLD

Transfer section		Cone
Retention in impression		Groove

Old cap for re-ordering:  
ø 3,3 mm, 5 pieces: 936030

### NEW

	Round with 2 faces
	Two wings

New cap:  
ø 3,3 mm, 5 pieces: 936034

### OLD

Transfer section		Cone
Gingiva/ Emergenzprofil		Round
Connection		HEX

Can no longer be ordered

### NEW

	Cone with locking function
	Round with 2 faces
	Prism-HEX

	GH 4.0 mm	GH 7.0 mm
4.2 mm	936026	936027

### OLD

Transfer section		Cone
Retention in impression		Groove

Old cap for re-ordering:  
ø 4,2 mm, 5 pieces: 936031

### NEW

	Cone with locking function
	Two wings

New cap:  
ø 4,2 mm, 5 pieces: 936035

As a service, we continue to offer the old transfer caps. However, we recommend that you replace the old transfer abutments and caps with new ones. Your area sales manager will be happy to offer you an attractive upgrade.